

MEDECO B.V.

Product Specifications

Brand : Klinion® Protection → ULTRA COMFORT

Product : Nitrile, powder-free, non-sterile, examination gloves

Packaging & Product Reference

Product Reference	Description	Size	Color	Box Dimension (mm)*	Piece(s) per Box	Outer Carton Dimension (mm)*	Box(es) per Outer Carton
102601	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	XS	White	235 * 110 * 75	150	390 * 230 * 245	10
102602	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	S	White	235 * 110 * 75	150	390 * 230 * 245	10
102603	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	M	White	235 * 110 * 75	150	390 * 230 * 245	10
102604	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	L	White	235 * 110 * 75	150	390 * 230 * 245	10
102605	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	XL	White	235 * 110 * 75	150	390 * 230 * 245	10
102606	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	XS	Blue	235 * 110 * 75	150	390 * 230 * 245	10
102607	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	S	Blue	235 * 110 * 75	150	390 * 230 * 245	10
102608	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	M	Blue	235 * 110 * 75	150	390 * 230 * 245	10
102609	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	L	Blue	235 * 110 * 75	150	390 * 230 * 245	10
102610	Klinion Personal Protection ULTRA COMFORT Nitrile Exam Glove	XL	Blue	235 * 110 * 75	150	390 * 230 * 245	10

*Sizes are approximate

Product description

Disposable, latex-free, powder-free examination gloves made of nitrile butadiene rubber.

Intended use

The Klinion Protection examination gloves are intended to be worn on the hand for medical purposes and to help prevent contamination between patient and examiner.

User Info available at www.medeco.org

Classification

MDD 93/42/EEC Class I

PPE Regulation (EU) 2016/425 Category III

→ EU declaration of conformity at www.medeco.org

MEDECO B.V.

Product Specifications

Features

- Economical and hygienic protection for hands
- Ambidextrous, suitable for left and right hand
- Finger-textured surface
- Powder-free

Product components

- Acrylonitrile Butadine Rubber (NBR)
- Products listed on this specification are:
 - Latex-free
 - Phthalate-free
 - Absent of human or animal tissue
 - Absent of blood or blood derivative

Specifications

Dimensions (mm)

Width	XS	75 ± 5
	S	85 ± 5
	M	95 ± 5
	L	106 ± 5
	XL	116 ± 5
Thickness	Cuff	Min. 0.04
	Palm	Min. 0.05
	Finger	Min. 0.06
Length	Standard Cuff	Min. 240

MEDECO B.V.

Product Specifications

Physical Properties

Unaged	Tensile Strength	Min. 14 MPa
	Elongation at Break	Min. 500 %
	Force at Break	Min. 6.0 N
Aged	Tensile Strength	Min. 14 MPa
	Elongation at Break	Min. 400 %
	Force at Break	Min. 6.0 N
Freedom from Holes	AQL 1.5	
Powder Residue	Max 2 mg/glove	

Quality Standards

- Manufactured in accordance with Quality Assurance ISO 13485 and ISO 9001

In compliance with:

- EN 420:2003 + A1:2009 (fit for special purpose)
- EN ISO 374-1:2016, EN 374-2:2015, EN 374-4:2013, EN ISO 374-5:2016
- EN 455-1:2000, EN 455-2:2015, EN455-3:2015, EN455-4:2009
- EN 16523-1:2015

Sterilization

N/A - Product is provided non-sterile

Storage conditions

Store in a cool and dry place. Opened boxes should be kept away from fluorescent and sunlight. Gloves are packed in dispenser which is suitable for transport. Keep the gloves in the box when not in use.

Shelf life

3 years

MEDECO B.V.

Product Specifications


Pictures



Ref. product 102602



Ref. product 102605

Approval	Name	Date	Signature
Product Manager	Levine van Leeuwen	20-12-2018	
RA-Manager	Francis Dekker	20-12-2018	