

Effective Pressure Area Care

MATTRESS OVERLAY

Material:

High grade, antimicrobial TPU material with optimal envelopment and immersion properties to ensure clinical efficacy and patient comfort.

Handwells:

The 'dog bone' handwells facilitate **easy log rolling**, mobilizing and patient transfers.

Indications for use:

- to aid in the prevention of pressure injuries
- can be used for all levels of risk
- to aid in the treatment of all categories of pressure ulcers (I - IV)

Easy Inflation:

- Clear valve for inflation and adjustments.
- Red CPR valve for **rapid deflation**.

Ventilation Holes:

- aiding the even distribution of pressure across the patient's body surface area.
- provides effective ventilation and optimises micro-climate.
- minimises risk of skin maceration as body fluids are drawn away from the patient's skin.

Security Straps:

Corner Straps secure and centre the Mattress Overlay on the bed.

4000A Mattress Overlay



4000 B includes the EHOB Double Action Air Delivery Pump

EHOB® products have been clinically evaluated as part of effective protocol for the prevention and treatment of pressure injuries.



EHOB® has provided effective solutions for prevention and treatment of pressure injuries for more than 30 years.

We offer an extensive range of highly effective pressure relieving products to fit your specific needs.

EHOB® Mattress Overlays examples of use:



Hospital



Homecare



Nursing Homes and Rehabilitation Centres

EHOB® Mattress Overlays:

- Highly cost effective
- Easy to use
- Latex free
- Low profile depth (2 to 5cm)
- One year warranty
- Up to 360 kg weight limit
- Can be inflated whilst the patient is lying on the mattress overlay
- Antimicrobial agent incorporated which is effective against a wide range of microorganisms including MRSA, C. Diff and E.Coli



4020A Theatre/Trolley Overlay



CL1200 CAIR Lite Pump

EHOB® Mattress Overlay Range

Order no	Size	Weight Limit	Description
4000A	193x86x8 cm	270 kg	TPU Mattress Overlay
4000B	193x86x8 cm	270 kg	As above with MAD pump
4010A	193x119x8 cm	360 kg	TPU Bariatric Overlay with MAD pump
4020A	198x56x5 cm	270 kg	Theatre/Trolley Overlay with MAD pump
546OR	Assure Surgery Table Straps 12/Bag (Straps inclusive for model 4020)		
540MAD	Handpump for Mattress Overlays & Riser Recliner Pad		
CL1200	CAIR Lite Pump for pressure checking and adjustments		

Can be reused following local policy guidelines for cleaning and disinfection. All EHOB products in the Static Air range incorporate an antimicrobial agent that acts as an effective barrier against the following microorganisms: Methicillin-Resistant Staphylococcus Aureus (MRSA); Clostridium difficile; Escherichia coli; Klebsiella pneumonia; Pseudomonas aeruginosa; Salmonella choleraesuis; Candida albicans; Enterococcus faecium; Aspergillus niger and Bacillus subtilis.



EHOB Europe - Lichtenauerlaan 102-120
3062 ME Rotterdam
T: +31 (0)10 7892404 - paul.mersel@ehob.com

ehob-europe.com